

COD. 6057/R – ENDOCAR
(endoscope trolley – dimensions 710 x 600 x 1350 h mm)



Body – Structure made of powder-coated steel RAL 7035, thickness 10/10, finished with acrylic polyurethane enamels resistant to scratches, stains and corrosion. Painting process complies with UNI EN ISO 9002 regulations. The structure can be washed and disinfected.

Worktop – Worktop made of self-extinguishing ABS, dimensions 520x400x60 h mm with raised edges on three sides. The worktop is resistant to scratches and stains and can be cleaned, washed and disinfected.

Drawer unit – Drawers are made of powder-coated steel and run on telescopic slides allowing their complete extraction. Drawers are provided with ergonomic and anti-trauma handles and stopper. Centralized key lock.

Available drawer fronts colours: red RAL 3000, green RAL 6029, sky blue RAL 5015, yellow RAL 1018.



Standard drawer configuration:

- 1 drawer, dimensions 425x400x100 h mm;
- 2 drawers, dimensions 425x400x150 h mm;
- 2 drawers, dimensions 425x400x200 h mm.

Left side – Ergonomic and anti-trauma ABS push handle

Right side – Cabinet made of powder-coated steel, thickness 8/10, dimensions 405x160x1220 h mm, equipped with 2 plastic supports for endoscopes. It is closed with 1 door with tempered glass, thickness 6 mm, provided with handle and key lock.

Base – Double perimetral and shock absorbent bumper to prevent the trolley from turning over. 4 thermoformed swivel castors, Ø 125 mm, mounted on bearings; 2 are lockable.