

Product specification

Bipolar coagulation forceps

Single Use, Disposable

Sterile disposable bipolar forceps with cable

- high gloss polished tip
- high-quality insulation
- Cable 3 m long
- Connector 2-pin international pin distance 22 to 28 mm



New

The insulation from the reusable product area; This blue insulation is better than the yellow one and is now also used for the disposable forceps.

The cable as well as the transition piece from tweezers to cable has also been improved and is made of black polymer.



Disposable bipolar coagulation forceps, „Jeweller“, straight, tip 0.5 mm, packing unit 12 pcs./box, sterile, single packed

Length
115 mm


Order-no.
241-SUBH-604805IC



Disposable bipolar coagulation forceps, „Adson“, straight, tip 1.0 mm, packing unit 12 pcs./box, sterile, single packed

Length
125 mm

Order-no.
241-SUBH-60601IC



Disposable bipolar coagulation forceps, straight, tip 1.0 mm, packing unit 12 pcs./box, sterile, single packed

Length
200 mm

Order-no.
241-SUBH-60571IC

Product specification

Bipolar coagulation forceps

Single Use, Disposable

Sterile disposable bipolar forceps with cable


- high gloss polished tip
- high-quality insulation
- Cable 3 m long
- Connector 2-pin international
pin distance 22 to 28 mm



New

The insulation from the reusable product area; This blue insulation is better than the yellow one and is now also used for the disposable forceps.

The cable as well as the transition piece from tweezers to cable has also been improved and is made of black polymer.




Disposable bipolar coagulation forceps, angled, tip 1.0 mm, packing unit 12 pcs./box, sterile, single packed

Length

180 mm

Order-no.

241-SUBH-60551AIC



Disposable bipolar coagulation forceps, bajonett, straight, tip 1.0 mm, packing unit 12 pcs./box, sterile, single packed

Length

200 mm

Order-no.

241-SUBH-60591IC